

Vedr.:  Mini Crosser batteribytte	PUP nr.:	24-007
	Udarbejdet af:	LPA
	Dato:	09-10-2024
	Intern <input checked="" type="checkbox"/>	Ekstern <input checked="" type="checkbox"/>

Vi har valgt at udskifte batteriet:

1523-1138 Batteri AGM 85-12 Power (MPX – MOVE)

Med det nye batteri:

1523-1157 Batteri MPX 90-12 (12V-100 Ah AGM)

De nye batterier er parret spændingsmæssigt, og er mærket med numre.

Det betyder at man ikke kan købe et ulige antal batterier. Altid kun lige antal.

Ændringen træder i kraft på nye køretøjer fra og med serienummer 03-57390-24-10.

Da mængden vi modtager lige nu, ikke kan dække behovet til både nye køretøjer og reservedelssalget, kan de nye batterier først købes som reservedele fra 6. januar 2025.



**Klaus Smedegaard**  
Administrerende direktør



Regarding:  Mini Crosser battery change	PUP No.:	24-007
	Prepared by:	LPA
	Date:	09-10-2024
	Internal <input checked="" type="checkbox"/>	External <input checked="" type="checkbox"/>

We have chosen to replace the battery:

1523-1138 Battery AGM 85-12 Power (MPX – MOVE)

With the new battery:

1523-1157 Battery MPX 90-12 (12V-100 Ah AGM)

The new batteries are matched according to voltage and are marked with numbers.

This means that you cannot buy an odd number of batteries. Always even numbers only.

The change comes into effect on new vehicles from and including serial number:

03-57390-24-10.

Since the amount we are receiving right now cannot cover the need for both new vehicles and spare parts sales, the new batteries can only be purchased as spare parts from 6 January 2025.



**Klaus Smedegaard**  
General Manager



medemagruppen